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Store between 15°C and 30°C (59°F and 86°F)
For *in vitro* use by professionals only.
We recommend reading entirely these instructions before using the test.

Introduction

This kit is a highly sensitive, specific and rapid lateral flow test that uses doxycycline-specific monoclonal antibody to detect its presence in urine samples. A positive test confirms the presence of doxycycline in the urine.

Basics of the test

This is a ready-to-use membrane test using red polymer particles labeling (test line) and blue polymer particles labeling (control line) that identifies the presence of doxycycline molecules in the urine sample.

If doxycycline is present in the urine sample at a concentration higher than the detection limit, the anti-doxycycline conjugates will react with these doxycycline molecules and **no test line will appear on the test strip.**

Urine will continue to migrate until it reaches a control reagent (upper line) that binds the control conjugate, thus producing a "migration control" line indicating that the test is functioning correctly.

Kit contents

- Disp&FLOW – Doxycycline kit contains the following:
 - 2 vials with desiccant cap each containing 25 strips for doxycycline testing
 - 50 disposable micropipettes (with a fixed volume of 400 µL)
 - 50 test tubes
 - Instructions for use.

Material necessary but not provided

- Gloves
- Sampling vials
- Timer (option)

Precaution

- All reagents are for *in vitro* detection use only.
- The kit can be used until expiry date under condition that it has been stored as indicated in this manual.
- Do not use the test after expiry date.
- Do not use the test if the packaging is damaged or if the vial is opened.
- The vials containing the strips must be kept at room temperature (between 15 and 30°C).
- All operations related to the use of this kit must be performed under strict hygienic conditions. Wear gloves when handling samples.
- Strips should be **handled only** by the upper orange part. Avoid any contact with the middle part of the strip which is composed of fragile nitrocellulose.
- The vials holding the strips contain a desiccant (in the cap) and should be kept tightly closed. The vial must be closed as soon as the required numbers of strips has been removed.** Strips are sensitive to ambient humidity – Do not use a strip more than 10 minutes after its removal from the vial.
- Use only tubes and disposable pipettes provided within the kit. Never use reagents from another kit (using other tube size and/or quality may affect the test).
- Do not dip the strips deeper in the urine sample than the black line under the arrows.**
- Strips and accessories provided within this kit are for single use only.**

Waste disposal

- Dispose of gloves, test tubes and used kits in accordance with biomedical waste regulations.
- Each user is responsible for the management of waste produced, and must ensure that it is disposed of in accordance with applicable legislation.

Storage and stability

- For optimal performance the kit should be stored between 15°C and 30°C (59°F and 86°F) in a dry place and protected from direct exposure to sunlight.
- Do not freeze the strips as this may cause malfunction of the test**
- The strips should be stored in provided desiccant tubes hermetically closed.**
- The kit expires on the date indicated on the box (or **6 months after opening** within the indicated expiry date)

Sample handling and collection

- Urine samples should be collected in a clean container.
- No centrifugation or filtration of urine is required.
- For optimal results urine samples must be tested as soon as possible after collection no more than 4 hours before being tested.**
- Alternatively, urine samples can be (i) stored between 2°C and 8°C (36°F to 46°F) before being tested or (ii) frozen and kept at -20°C (-4°F) for later use. In all cases, samples must be brought to room temperature before testing.
- Disperse any precipitate by gentle agitation (prolonged storage).
- Avoid excess exposition to light as doxycycline may be degraded** (cover with aluminum paper for optimal protection if necessary).
- Urine samples should be handled with care as they are potentially infectious.

Procedure (cf. figure 1)

Before testing, allow the urine samples to reach room temperature (15°C – 30°C or 59°F – 86°F) in their closed container.

- Place 400 µL (+/-50) of urine sample in the test tube
 - Use a disposable pipette
 - Or use a laboratory pipette
- Remove a strip from the desiccant vial.
- Place the strip vertically into the test tube containing the urine sample (**do not immerse the strip above the line placed under the arrows**).
- Incubate for 8 +/- 2 minutes.
- Remove test strip from the urine sample and immediately read the results.

Ignore the appearance of any new line appearing on the test strip after the waiting time.

Reading should be done on humid strip after 10 minutes maximum.

Interpreting results

The results are to be interpreted as follow (cf. Figure 1):

1- **Negative test = No doxycycline has been detected:** Two (2) lines are present on the test strip. When the upper control line (blue) and the lower test line (red) appear on the test strip, the test must be considered as **NEGATIVE. The urine sample does not contain doxycycline.** Any test line (red), even weak, should be considered as a negative result.

2- **Positive test = Doxycycline has been detected:** Only one (1) line is present (blue control line). The test must be considered as **POSITIVE. The urine sample contains doxycycline.**

3- **Invalid test:** There is no detectable upper blue control line. The absence of a control line indicates a failure in the test procedure. **A new test must be performed** using a new urine sample and new test strips.

Note: After the drying process, a very faint shadow may appear at the test line position. It should not be regarded as a negative result.

To store the result, allow the test strip to dry after removing the absorbent material at its base.

Performance

A. Detection limit

The detection limit of the test is $\geq 5\mu\text{g/mL}$ Doxycycline. The test provides POSITIVE results (no test line, only the upper blue control line) for doxycycline concentrations above the detection limit.

B. Evaluations

All field tests demonstrated a positive result with all the strips tested (presence of doxycycline) when the doxycycline had been taken. The presence of doxycycline in urine samples was confirmed in a specialized lab by HPLC (High Performance Liquid Chromatography)

C. Accuracy

In order to check intra-batch accuracy, the same positive and negative samples were processed 20 times on test strips of the same production batch under the same experimental conditions. Results were confirmed for 100% of the samples.

In order to check intra-batch accuracy, some samples (positive and negative) were processed on the strips from three (3) different production batches. All results were confirmed to be accurate.

D. Interference

The test shows interference with the following antibiotics: Anhydrotetracycline, Méthacycline (LOD<5µg/ml), 4-Epitetracycline (LOD<50µg/ml) and Oxytétracycline (LOD<100µg/ml). But all laboratory tests demonstrated no cross-reactivity with respect to the following molecules at 100µg/mL: Proguanil, Mefloquine, Minocycline, Chlorotetracycline and Tetracycline.

Limitations of the kit

This test is only qualitative and cannot measure accurately the concentration of Doxycycline in urine samples.

The test strips cannot determine whether effective protection levels against malaria have been obtained.


They can only confirm antimalarial drug intake within the limit of detection.

Technical problems / Complaints

If you encounter a technical problem or if results do not correspond to those indicated in the instructions:

- Note the batch number of the incriminated kit.
- If possible, store the problematic test strip at room temperature, and urine sample at -20°C (-4°F)
- Contact BIOTEM (info@biotem.fr) or your local distributor

Figure 1

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