

BIOTEM  
885, rue Alphonse Gourju  
38140 Apprieu  
FRANCE

Tel: +33 (0)4 76 65 10 91

Fax: +33 (0)4 76 67 48 86

[info@biotem.fr](mailto:info@biotem.fr)

[www.biotem-antibody.com](http://www.biotem-antibody.com)

.....TECHNICAL DATASHEET.....

## Human IgG4 Isotype Control Ref # BIO.061.0

### Product description

Description	Monoclonal antibody to a synthetic hapten (polychlorobiphenyl). This antibody is a negative control.
Isotype	Human IgG4, kappa recombinant antibody
Specificity	The antibody is directed against a chemical small molecule. Do not cross react against mammalian protein.

### Properties

Form	Liquid
Storage instructions	Store at +4°C short term (< 1 month). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles
Storage buffer	Phosphate Buffer 10mM - NaCl 0.15M - pH 7,4
Production	Production in CHO cells
Purification	Purified IgG fraction prepared by affinity chromatography on protein A
Endotoxin Level	<10 EU/mg by the LAL method (see certificate of analysis for the level of each lot)
Restrictions	For research use only. Not for use in diagnostic or therapeutics procedures.

### Applications

ELISA	Use at an assay dependent dilution
In vivo	Use at an assay dependent dilution
Flow cytometry	Use at an assay dependent dilution
Optimal dilutions/concentrations should be determined by the end user	