

BIOTEM
885, rue Alphonse Gourju
38140 Apprieu
FRANCE

Tel: +33 (0)4 76 65 10 91

Fax: +33 (0)4 76 67 48 86

info@biotem.fr

www.biotem-antibody.com

.....TECHNICAL DATASHEET.

Human IgG1 silent Isotype Control Ref # BIO.063.0

Product description

Description	Monoclonal antibody to a synthetic hapten (polychlorobiphenyl). This antibody is a negative control.
Isotype	human IgG1 LALA ASPS, kappa recombinant antibody
Specificity	The antibody is directed against a chemical small molecule. Do not cross react against mammalian protein.

Properties

Form	Liquid
Storage instructions	Store at +4°C short term (< 1 month). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles
Storage buffer	Phosphate Buffer 10mM - NaCl 0.15M - pH 7,4
Production	Production in CHO cells
Purification	Purified IgG fraction prepared by affinity chromatography on protein A
Endotoxin Level	<10 EU/mg by the LAL method (see certificate of analysis for the level of each lot)
Restrictions	For research use only. Not for use in diagnostic or therapeutics procedures.

Applications

ELISA	Use at an assay dependent dilution
In vivo	Use at an assay dependent dilution
Flow cytometry	Use at an assay dependent dilution
Optimal dilutions/concentrations should be determined by the end user	