

**BIOTEM Apprieu – Headquarters**  
**Hybridoma & Recombinant**  
**Antibody Departments**  
Parc d'activités Bièvre Dauphine  
885, rue Alphonse Gourju  
38140 Apprieu  
France

**BIOTEM Colombe**  
**Immunoassay Department**  
Parc d'activités Bièvre Dauphine  
100, rue Eugène Perrin  
38690 Colombe  
France

**Contact**

Commercial: +33 (0)4 76 93 97 35  
General: +33 (0)4 76 65 10 91  
Email: [info@biotem.fr](mailto:info@biotem.fr)

[www.biotem-antibody.com](http://www.biotem-antibody.com)

**TECHNICAL DATASHEET**

# Human IgG1 silent Isotype Control

## Ref # BIO.063.0

## Product Description

Description	Recombinant antibody to a synthetic hapten This antibody is a negative control.
Isotype	Human IgG1 LALA ASPS, kappa recombinant antibody
Specificity	The antibody is directed against a chemical small molecule. Do not cross react against mammalian proteins.

## Properties

Concentration	Refer to the certificate of analysis
Storage instructions	Store at +4°C short term (1-2 month). Aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles
Storage buffer	Phosphate Buffer 10mM - NaCl 0.15M - pH 7.4
Production	Production in CHO cells
Purification	Purified IgG fraction prepared by affinity chromatography on protein A
Endotoxin Level	<10 EU/mg by the LAL method (see certificate of analysis for the level of each lot)
Restrictions	For research use only. Not for use in diagnostic or therapeutics procedures.

## Applications

ELISA	Use at an assay dependent dilution
In vivo	Use at an assay dependent dilution
Flow cytometry	Use at an assay dependent dilution
Optimal dilutions/concentrations should be determined by the end user	