

Disp&FLOW[®] FOOD

Rapid Test for Food Industry



Disp&FLOW[®]

Rapid tests for the detection of allergens or prohibited components on your production lines






Overview

Our **Disp&FLOW®-Food** product is a rapid immuno-chromatographic test for the qualitative determination of different components in food, or on cooking utensils and preparatory surfaces in the context of **food production** and **relevant processing facilities**.

This product line is therefore aimed at food industry professionals for first-line rapid control tests.



References available

Food Allergen / Intolerance - Prohibited Components

Reference No.	Name	Reference No.	Name
BIO.068	Sesame	BIO.035	Almond
BIO.067	Walnut	BIO.034	Peanut
BIO.049	Frog	BIO.033	Hazelnut
BIO.048	Macadamia Nuts	BIO.032	Fish
BIO.047	Cashew Nuts	BIO.031	Milk
BIO.046	Brazil Nuts	BIO.030	Egg (Ovotransferrine)
BIO.045	Celery	BIO.029	Egg (Ovalbumin)
BIO.040	Gluten	BIO.028	Shrimp
BIO.039	Garlic	BIO.027	Pork (Ready)
BIO.038	Mustard	BIO.026	Pork (Blood / fat)
BIO.037	Soy	BIO.023	Human Blood
BIO.036	Coconut	BIO.022	Alcohol

*Ask our team
for other references*



BIOTEM is specialized in the generation and production of **custom antibodies and immunoassays** (ELISA & Lateral Flow rapid tests) since 1980. Our team can assist you in the **development of your own product** according to your specifications or can expand its product range according to the demand.

Contact us to discuss your needs! >> INFO@BIOTEM.FR